# **GBTx** Communication Document

University of Maryland LHCb group

April 19, 2018

Rev: 0.3 (cf53be3)

### Contents

1	Har	dware setup	2				
	1.1	Overview	2				
	1.2	Configure GBTx to use external $\mathrm{I}^2\mathrm{C}$ adapter	2				
	1.3	Configure slave GBTx to use master SCA channel $\dots$	2				
	1.4	Configure GBTx to use GBT-IC channel	3				
	1.5	Reset GBTx	4				
			4				
2	Software setup						
	2.1	Program master GBTx with external $\mathrm{I}^2\mathrm{C}$ adapter	4				
	2.2	Check communication between master GBTx and MiniDAQ $$ .	5				
	2.3	Program individual registers of slave GBTx	6				
	2.4	Program slave GBTx with configuration files	7				
$\mathbf{A}$	ppen	dices	8				
	A	Fix "Waiting for a GBT server to run"	8				
	В	Bypass components in GBT data path	8				

### 1 Hardware setup

#### 1.1 Overview

Our current setup consists of one master and one slave GBTx board. The master is connected to the MiniDAQ GBTx channel 3 (fiber 8), and the slave can be connected to either GBT channel 0 (fiber 6) or channel 6 (fiber 11). The master synchronizes its on-board clock to the signal from the MiniDAQ, and propagates its clock signal to the slave. The slave does not have an on-board clock, and is configured to obtain clock signal externally.

The master I<sup>2</sup>C port is connected to an external USB device. The slave I<sup>2</sup>C port is connected to the master. Both are set to be programmed by the I<sup>2</sup>C channel, rather than GBT-IC channel.

The current setup is capable of:

- 1. Program the slave GBTx board with MiniDAQ directly.
- Read/Write the register value of the master GBTx board with GBT-IC specification on the MiniDAQ.
- 3. Do PRBS tests from MiniDAQ to the slave, then back to the MiniDAQ. The master is also required as the slave can only obtain its reference clock from the master.

### 1.2 Configure GBTx to use external I<sup>2</sup>C adapter

This setup is required to program a GBTx board using an external I<sup>2</sup>C adapter. Follow Figure 1 to connect an external I<sup>2</sup>C adapter.




Figure 1: Schematic for external I<sup>2</sup>C adapter setup.

The  $I^2C$  adapter used in our lab is made in-house. For the 2x1 connector, make sure the side that has *no* metal contact is facing up.

#### 1.3 Configure slave GBTx to use master SCA channel

This setup is required to program a slave via a master SCA channel. In a typical scenario, the master is connected to a MiniDAQ so that programming the slave using the MiniDAQ directly<sup>2</sup> is possible. Follow Figure 2 to connect a slave GBTx to the SCA channel of a master GBTx board.

<sup>&</sup>lt;sup>1</sup> The master GBTx is also connected to the MiniDAQ with a different channel, to provide reference clock to the slave.

 $<sup>^2~{\</sup>rm MiniDAQ} \rightarrow {\rm master~GBTx} \rightarrow {\rm slave~GBTx}.$ 




Figure 2: Schematic for slave to master SCA setup.

The black ground cable can be connected to any of the ground pin on the master GBTx.

There is a 2x1 to 2x1 cross-type cable made in-house to replace the redblue cables. To use that cable, make sure the two 2x1 connectors have the same orientation (e.g. the sides *without* metal contact are both facing up).

### 1.4 Configure GBTx to use GBT-IC channel

It might be useful to read/write individual registers from/to a GBTx board. In this case, follow the Figure 3 to flip the configSelect switch.

Flip the configSelect switch will render the external I<sup>2</sup>C adapter ineffective. None of the GBTx register value is fused onto the board, so a GBTx board in our lab must always be programmed externally via I<sup>2</sup>C before flipping the switch.




Figure 3: Schematic for flipping the configSelect switch. A jumper should be used to connect the two pins marked above.

#### 1.5 Reset GBTx

Sometimes GBTx boards will not be properly reset by reprogramming. In such case, a hard reset is needed. Follow the Figure 4 to reset GBTx boards.



Figure 4: Schematic for resetting GBTx boards. A jumper should be used to connect the two pins marked above.

## 2 Software setup

### 2.1 Program master GBTx with external I<sup>2</sup>C adapter

Before proceed, follow the instruction on subsection 1.2 to configure the hardware first. As said in the previous section, the master GBTx board must be programmed via an external  $I^2C$  adapter first. A Windows 7 computer on the rack is used. The **GBTX programmer** is located at:

DT\_Rack\GBTx\_programmer\GBTxProgrammer.jar

SB dongle open	ed succe	essfully - f	w v1.c	Register S	Select	#Register	Value (hex)	Name	
15	)			0	0x00	ckCtr0			
1 sent to VTarget LDO GBTX found on address 1						1	0x00	ckCtr1	
						2	0x00	ckCtr2	
			3	0x00	ckCtr3				
Selected GBTX address 1						4	0x00	ttcCtr0[7:0]	
						5	0x00	ttcCtr1[7:0]	
Configuration loaded from txt file						6	0x00	ttcCtr2[7:0]	
•				7	0x00	ttcCtr3[7:0]			
ritton and road	rogictor	250 00113	I Programming			8		ttcCtr4[7:0]	
Written and read registers are equal. Programming						9	0x00	ttcCtr5[7:0]	
as successful!				10	0x00	ttcCtr6[7:0]			
						11	0x00	ttcCtr7[7:0]	
						12	0x00	ttcCtr8[7:0]	
						13	0x00	ttcCtr9[7:0]	
						14 15	0x00 0x00	ttcCtr10[7:0]	
								ttcCtr11[7:0]	
				<b>V</b>		16	0x03	ttcCtr12[7:0]	
	<b>V</b>		17	0x03	ttcCtr13[7:0]				
Load GBTX configuration				V		18	0x03	ttcCtr14[7:0]	
				<b>V</b>		19	0x03	ttcCtr15[7:0]	
Write ALL to the GBTX				V		20	0x03	ttcCtr16[7:0]	
				<b>V</b>		21	0x03	ttcCtr17[7:0]	
				V		22	0x03	ttcCtr18[7:0]	
Write selected to the GBTX						23	0x03	ttcCtr19[7:0]	
						24	0xff	ttcCtr20[7:0]	
ъ.			_	~		25	0x03	ttcCtr21[7:0]	
Read all registers					V	26	0x7f	ttcCtr22[7:0]	
				V		27	0x28	serCtr0[7:0]	
Reset 1V	5	Re	set 2V5			28	0x00	txCtr0	
				V		29	0x15	txCtr1	
BTX on 1				V		30	0x15	txCtr2	
	Rea	d state:	24 (dec)	V		31	0x15	txCtr3	
	Statue:	ldle (norm	al etatue	V		32	0x66	txCtr4	
	Status: Idle (normal status when running)			V		33	0x00	txCtr5	
	wnen ri	when running)				34	0x0d	desCtr0	
	- Down	Power GBTX trough I2C a		V		35	0x42	rxCtr0	
	Powe	I GD I X LIOU	yıı izc avapter			36	0x00	rxCtr1	
	E	nable exp	ert mode	V		37	0x0f	rxCtr2	
				V		38	0x04	rxCtr3	
		V		39	0x08	rxCtr4			
Fuse selected registers					ect all	registers		DEselect all registers	

Figure 5: A typical UI for GBTX programmer.

Launch the programmer, a typical UI is shown in Figure 5, Click Load GBTX configuration and load a configuration file, which is located at:

```
{\tt DT\_Rack\backslash GBTx\_programmer\backslash GBTx\_TRx\_v12\_test\_with WatchDog.txt}
```

Then click Write ALL to the GBTX. Check the returned message to make sure everything works (supposedly). Now click Read state. If the master GBTx is configured correctly and is connected to a working MiniDAQ, the return value should be:

```
24 (dec): Idle (normal status when running)
```

#### 2.2 Check communication between master GBTx and MiniDAQ

After program the master GBTx board and verify the return value, we can check the communication between the GBTx and the MiniDAQ with **GBT** Client.




Figure 6: A typical UI for GBT Client, SCA tab.

Here we need to use the Linux server on the rack. To launch that program, locate the **gedi** panel. Under **JCOP** framework, click the **GBT** Client. Choose **GBT** option under the **Communication** tab. Navigate to **SCA** tab.

To check whether the link between GBTx and MiniDAQ is successfully established, configure the parameters *exactly* as shown in Figure 6. Now click **Write/Read**. The **Data out** field should have a value of "FE", and the **Status** should be "OK (Write/Read)".

**GBT ID** corresponds to the physical optical link that is connected to the master GBTx board. Recall that in our setup, the master is connected to optical fiber 8, which is mapped to GBT channel 3.

There are only 4 **Registers**. All 4 registers work, the "Control Register B" is chosen for no apparent reason. What matters is the return value should be consistent with our expectation.

### 2.3 Program individual registers of slave GBTx

Before proceed, follow subsection 1.3 to set up hardware. Again we use **GBT Client** to program individual registers of the slave GBTx. Configure the parameter *exactly* as shown in Figure 7.




Figure 7: Correct parameters to configure GBT I<sup>2</sup>C tab.

If the Write command failed, try Activate Channel first.

**Address** corresponds to GBTx index, which is configured via a onboard DIP switch. Default to 1.

Recall that 1 Byte is 8 bits, and 1 Byte corresponds to a 2-digit hex number.

Notice in the **Data out** field, the last 3 numbers are the same. This is a requirement in the GBT specification that the configuration registers must be triplicated to minimize error.

### 2.4 Program slave GBTx with configuration files

The setup is exactly the same as shown in subsection 2.3. All we need to change are replace the **Size** field with "0", and provide the configuration file location with the following syntax, which needs to be filled in the **Data** in field:

file: <path\_to\_configuration\_file>

This is equivalent to program the slave with an external I<sup>2</sup>C adapter. The advantage of doing this way is we don't need to remove the I<sup>2</sup>C adapter from the master. In other words, once the hardware setup is done, we don't have to touch them anymore; whereas if we program both master and slave with an I<sup>2</sup>C adapter, it must be moved around between the master and the slave.

### **Appendices**

### A Fix "Waiting for a GBT server to run"

Sometimes when we launch the **GBT Client**, it would claim that it is "Waiting for a GBT server to run". Take this warning literally: It means that currently there is no GBT server that is alive. The fix is easy: fire up a terminal, and type in "GbtServ".

### B Bypass components in GBT data path

As shown in Figure 8, we can bypass components in GBT input/output path.




Figure 8: GBT data path. Stolen from LHCb GBTX manual.

The register 0x1c and 3 subsequent registers control the output path; 0x2d controls the input path. Each of the paths have 3 configurable selectors, each has 4 operations modes  $(0_{10} - 3_{10}^3)$ , which translate to  $0_2 - 11_2$ ). According to Figure 8, the control sequence should be given in:

Here is an example. Let us consider the Tx normal operating case (i.e. all 3 switches are configured to 1). We have:

$$\begin{array}{cccc} \underbrace{1_{10}}_{\text{TXSwitch2}} & \underbrace{1_{10}}_{\text{TXSwitch1}} & \underbrace{1_{10}}_{\text{TXSwitch0}} \\ & & & & \\ \underbrace{01_2}_{\text{TXSwitch2}} & \underbrace{01_2}_{\text{TXSwitch1}} & \underbrace{01_2}_{\text{TXSwitch0}} \end{array}$$

To convert binary  $010101_2$  to hexadecimal, we left pad the result with 2 additional 0's:  $0001,0101_2$ . The end result is  $15_{16}$ .

<sup>&</sup>lt;sup>3</sup> The subscript indicates the base.

 $15_{16}$  is exactly the last 3 bytes in  $\bf Data$  out field in Figure 7. This is not a coincidence.

Recall that each 4 digits in binary exactly correspond to 1 digit in hexadecimal.